

New Patient Details Feature on ZEPHYRx Provider Dashboard

In the latest update to the Provider Dashboard, ZEPHYRx has added a 'Patient Details' tab to the patient profile. Patient Details has consolidated existing Dashboard features while adding new functionality. The Patient Details tab is now where Providers will set or edit patient demographics, update a patient study ID (if applicable), assign a primary physician and enter patient notes. ZEPHYRx has also released a reminders and tasks feature that will allow Providers to assign tasks and set reminders for their patients through the Provider Dashboard. From the Patient Details tab Providers can send their patients reminders to complete FVC/SVC tests or games and monitor the status of their tasks in the Dashboard. Reminders/tasks is disabled by default but can be enabled at your organization's request.

[← Back to Patient List](#)

Patient Details

Breathing Exercises

Pulmonary Function Tests

Start Remote Session

Patient Details

+ Edit

Identification

ID2da853f3-ead2-46e8-b64d-4e2e8eb54f63

StudyIDSerkules13I

CFF ID9744860

Spirometer IDSM-005-Z068447

Demographics

Sex at BirthMale

Birth1988, 33 years

Height73.00 in

Weight225.00 lbs

Operational Settings

Ignored In OverreadNo

Primary Physician in ClinicNot Selected

Reminder Settings

FVC remindersMonthly

SVC remindersDaily

Breathing Exercise remindersDaily

Notes

Task Details

+ PFTnow!

Type	Done	Created	Expires	Remove
Game	✗	09/09/2021	09/10/2021	--
Pft: SVC	✗	09/09/2021	09/10/2021	--
Pft: FVC	✓	09/01/2021	10/01/2021	--

Edit Settings for Patient

Demographics

Name:

Sex at Birth:☒ Male ☐ Female

Birthdate:

6 ▾ / 2 ▾ / 1962 ▾

Ethnicity:South East A ▾

Height:67.00 in.

Weight:136.00 lbs.

Edit for selected date range:

Start Date:09/15/2021

End Date:09/15/2021

Operational Settings

Ignored In Overread:☐ Yes ☐ No

Select Primary Physician:Not Selected ▾

Reminder Settings

Send FVC Reminders:Never ▾

Send SVC Reminders:Never ▾

Send Exercise Reminders:Never ▾

Caution: Editing demographics will cause a recalculation of the patient's predicted values within the dates selected. This process may take several minutes.

If this patient has their Breathe Easy app open, please have them close and quit the app when changing demographics to see the changes.

Save changes

Cancel

If you would like to learn more about the Patient Details section or other Dashboard updates please message us at support@zephyrx.com or call our live support line at 888-452-6269. Thank you!